



Epithelix

in vitro Solutions for Respiratory Diseases and Chemical Testing




Streptococcus pneumoniae inhibits *Pseudomonas aeruginosa* growth on nasal human epithelium *in vitro*

Song Huang, Carole Bertinetti, Ophélie Verbeke, Mireille Caul-Futy, Paul Alouani, Ludovic Wiszniewski, Samuel Constant.
Epithelix, 18 Chemin des Aulx, Plan-Les-Ouates, CH-1228, Geneva, Switzerland

Pathogens colonizing the respiratory tract compete with a range of other bacteria. *Pseudomonas aeruginosa* (PA) infections are increasingly associated with acute exacerbations in chronic obstructive pulmonary disease. *Streptococcus pneumoniae* (SP), meanwhile is a main cause of pneumonia, meningitis, it can lead to infections and other respiratory diseases such as bronchitis. We report herein the use of 3D airway epithelia reconstituted *in vitro* to study interactions of PA and SP on nasal mucosa. MucilAir™, a fully differentiated human airway epithelium made of a mixture of primary nasal cells from 14 donors, was used to study the effects and behaviour of PA and SP (inoculated at $3E^{+02}$ and $3E^{+11}$ CFU/cm² respectively) cultured separately or together over 24 hours.



Testing Strategy



MucilAir™, a fully differentiated 3D *in vitro* cell model of the human airway epithelia. The cells were freshly isolated from the biopsies. After 45 days of culture at air-liquid interface, the epithelia were fully differentiated, both morphologically and functionally.

Information from Apical Side


- TEER measurement
- Enumeration of colony (CFU)
- Cilia beating Frequency monitoring
- Defensin / Cathelicidin & Mucin Secretion




Pseudomonas Aeruginosa (PA) and *Streptococcus Pneumoniae* (SP) were inoculated respectively at $3E^{+02}$ and $3E^{+11}$ CFU/cm² on MucilAir™-pool of donors on apical side. Toxicity induced by infection is monitored by several non-destructive endpoints such as TEER measurement, cilia-beating frequency monitoring, mucin secretion, IL-8 and LDH releases

General protocol and schedule of end point measurement, N= 3 per condition and per time point. Statistical comparison was done by Student's t test (*p<0.05).

Microscopic analysis



Scanning electron microscopy image of *Pseudomonas aeruginosa* (blue arrow) and *Streptococcus pneumoniae* (orange arrow) on MucilAir™ apical surface at 24 hours post inoculation during a concomitant infection.




Transmission electron microscopy images representing A) internalized bacteria of PA and SP on the apical surface, B) SP intracellular and C) migration of PA through the pores of insert membrane

Apical, basolateral and intratissular PA and SP growth were quantified by Colony Forming Unit (CFU). Impairment of epithelial homeostatic barrier function was evaluated through monitoring the Trans Epithelial Electrical Resistance (TEER), cytotoxicity (LDH), cilia activity, mucin and IL-8 release.


Results

Evaluation of barrier function and cytotoxicity on MucilAir™ after PA or SP infection alone or together after 24h




SP infection had no impact on TEER or LDH release, whereas PA induced a decrease of TEER coupled with cytotoxicity. In the co-infection model, only a drop of TEER is detected.

Quantification of PA and SP growth alone or together on MucilAir™ after 24h



After 24h of co-infection, inhibition of PA growth was observed on the apical surface

Soluble factors release on MucilAir™ after PA or SP infection alone or together after 24h



SP infection induced a strong mucin and H₂O₂ secretion whereas PA infection induced secretion of β-Defensin 2, Cathelicidin and IL-8.

PA infection induces a loss of TEER, 20% cytotoxicity, an increase of IL-8, and an up regulation of β2-defensin and cathelicidin. On the contrary, SP strongly increases the mucin production and H₂O₂ release. While inoculated together, a lower apical PA growth is observed ($- 3E^{+3}$ CFU/cm²) suggesting an inhibition due to the presence of SP.

Conclusion

These results suggest that *in vitro* human airway epithelia is a useful model to study bacterial interaction on the human nasal mucosa.

